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Kintex Kinesiologie Tape

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EN INSTRUCTIONS FOR USE

Kintex Kinesiologie Tape

Dear customer,
We are delighted that you have chosen Kintex Kinesiologie Tapes Classic. If you have
any questions, please contact our Customer Service via our website: www.kintex.de

For use as pain, muscle, lymph and correction tape.
Thanks to the wavy structure of its hypoallergenic acrylic adhesive, the skin-friendly
Kintex Tape is very adaptable and stretchable (40 - 60%). Our tape guarantees
and enables perfect adjustment to structures of the human body such as the skin, the
fascia, the muscles, the meridians and more. Its texture is very durable, waterproof and
tear-resistant.

Explanation of symbols

- Latex-free
- Do not use twice, adhesive strength decreases
- Protect from sunlight
- Read and follow the instructions for use
- Keep dry
- Medical device
- Please dispose of properly

Application

- Measure the desired length and cut the tape with scissors.
Please ensure that the corners are rounded.
- Remove the tape from the backing paper without reaching into the adhesive surface by
sticking on the base (approx. 5cm) and then sticking on the rest of the tape.
- Before applying the tape, the body is pre-stretched, this varies depending on the goal
and type of application.
- You can use the Kintex Tape for sports injuries, joint pain, back pain, oedema,
muscle tension and other indications and contraindications.*
- We offer you over 40 online instruction videos for various application techniques.

For more information visit our
website or scan
the QR code:

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Warning

Never use the tape on open wounds or as bandages. If you feel itching or an
allergic reaction, remove the tape. Before each application, a detailed diagnosis forms
the basis of a targeted therapy. We recommend having our tape applied exclusively
by trained users or doctors. *The effect of kinesiology tape has not been scientifically
proven, our statements are based on our many years of experience.

Material

97% cotton, 3% elastane, hypoallergenic acrylic adhesive

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